



Sure Quality is on  
*Sure-endo*


Taper .02




Taper .03




Taper .04




Taper .05




Taper .06




Protaper Size




Waveone Gold Size




Reciproc Size




Accessory Size (Regular)



Japanese Size



*Gutta Percha Points*



## *Gutta Percha Points* *Absorbent Paper Points*

Hand-Rolled Gutta Percha Points and  
Absorbent Paper Points for Root Canal from  
**Sure-endo**



*Just try it*  
feel the differences

Clean | Size | Variety  
Sterilized | Size | Variety




### Sure Dent Corporation

#809, 52, Sagimakgol-ro, Jungwon-gu, Seongnam-si, Gyeonggi-do, Rep. of Korea  
Tel. +82-31-737-4720 E-mail. [suredent@suredent.com](mailto:suredent@suredent.com) Web. [www.suredent.com](http://www.suredent.com)



SD-01 | Rev. 10 / 2019



# Gutta Percha Points

Hand-Rolled Gutta Percha Points

## Clean

Neat cutting

## Size

Accurate size of tip & taper

## Variety

Various sizes and packages



## Regular GP Points



## Bio GP Points



## Bioceramic Impregnated and Coated GP Points

Regular GP points in various sizes and packages of your choice

### Advantages

- Accurate size : Strict standard than ISO for top quality  
: Consistent quality  
: Laser inspection for tip & taper
- Vivid length mark : Calibrated at 16mm, 18mm, 19mm, 20mm, 22mm, and 24mm
- Great radiopacity
- Various size selections in diverse packages

Innovative Bioceramic Impregnated GP Points

### Advantages

- Biocompatibility
- Higher radiopacity than regular GP : 30-40% higher
- Clean and neat cutting with a trimmer : Not Sticky at all
- Stiff but flexible at the same time
- Great to use with sealers of any kinds : Intense bonding action
- Vivid length mark : Calibrated at 16mm, 18mm, 19mm, 20mm, 22mm, and 24mm

GP which is composed of bioceramic materials, creates intense-bonding obturation in canal wall

### Advantages

- The best use with bioceramic sealer : Coated with bioceramic material  
: The most powerful bonding forms when it used with bioceramic sealer
- Biocompatibility
- Stiff but flexible at the same time
- Higher radiopacity than regular GP : 30-40% higher
- Clean and neat cutting with a trimmer : Not Sticky at all

# Absorbent Paper Points

Sterilized endodontic Paper Points which is used to dry the canal and introduce medicaments

## Sterilization

Sterilized by Gamma-ray

## Size

Accurate size of tip & taper

## Variety


Various sizes and packages



## Paper Points

### Advantages

- High solidity and absorbency for operators' convenience
- 100% Irradiated by gamma-ray : Sterilization assurance level: 10<sup>-6</sup>
- Accurate size : Consistent quality  
: Laser inspection for tip & taper
- Lintless clean surface
- Various size selections in diverse packages
- Vivid length mark : Calibrated at 16mm, 18mm, 19mm, 20mm, 22mm, and 24mm



## Obturator

GP Obturator bullet for obturation gun

### Advantages

- Specially formulated
- Easy to heat and handle
- Biocompatibility
- Various types : Soft / Regular / Hard
- Highly flowable

### Possible Types \* Melting point

- Soft type : 150℃
- Regular type : 200℃
- Hard type : 230℃

### Package

- 100pieces

## Size Chart

Size	Taper	ISO Standardized Size				
		Taper .02	Taper .03	Taper .04	Taper .05	Taper .06
● 08		v				
● 10		v				
○ 15		v	v	v	v	v
● 20		v	v	v	v	v
● 25		v	v	v	v	v
● 30		v	v	v	v	v
● 35		v	v	v	v	v
● 40		v	v	v	v	v
○ 45		v	v	v	v	v
● 50		v	v	v	v	v
● 55		v	v	v	v	v
● 60		v	v	v	v	v
● 70		v	v	v	v	v
● 80		v	v	v	v	v
○ 90		v				
● 100		v				
● 110		v				
● 120		v				
● 130		v				
● 140		v				
Assorted		15 - 40	20 - 45	20 - 45	20 - 45	20 - 45
		45 - 80				
		90 - 140				

※ Both Length marked & Non-length marked are possible  
※ GP & PP are available in all of the above sizes

Type	Non-Standardized Size			
	Protaper	Protaper Next	Waveone Gold	Reciproc
Size	● F1	● X2	● Small	● R25
	● F2	● X3	● Primary	● R40
	● F3	● X4	● Medium	● R50
	● F4	● X5	○ Large	
	● F5			
	F1 - F3 F4 - F5	X2 - X5	S - L	R25 - R50

※ Both Length marked & Non-length marked are possible  
※ GP & PP are available in all of the above sizes

Type	Non-Standardized Size	
	Japanese	De-tray
Size	S	XXF
	FS	XF
	FM	F
	M	M
	FL	
	L	
	S - FL	

※ Only GP correspond to all of the above sizes

Type	Non-Standardized Size	
	Regular Type	
Size	GP	PP
	XF	XXF
	FF	XF
	MF	F
	F	M
	FM	C
	M	XC
	ML	XXF - XC
	L	
	XL	
	MF - M	

## Size Color Type

15	45	90	Large			
20	50	100	Small	F1 / F5	X5	R50
25	55	110	Primary	F2	X2	R25
30	60	120	08	F3	X3	
35	70	130	Medium			
40	80	140	10	F4	X4	R40

## Package Style



Slide Box for GP/PP



Round Vial for GP/PP



Endo Box for GP/PP



4D Box for GP



Square Vial for GP/PP



Cell Pack for PP

Sure-endo manufactures all sizes of GP & PP and meets the needs of all customers